



Stony Brook University

Surface Roughness Analysis for Atomic Force Microscope using Python

CME 502 Final Project Presentation



5/4/2023

Huiting Luo



Introduction

- Atomic force microscopy (AFM)
 - AFM is a widely used surface analysis technique for investigating nanoscale structure of materials.
 - Our lab has utilized a Bruker AFM for over 40 years and has been used in numerous research projects.
 - Explore the growth of cells on polymers scaffold
 - Examine the effects of polymer additives on perovskite grain boundary theory
 - Investigate the interaction between polymer thin films at the interface
- Surface roughness is an important parameter in characterizing the texture and visual qualities of a material's surface





Problem

- Bruker's NanoScope Analysis is a powerful software for processing AFM data
 - But unfortunately, it was only installed on certain outdated computers in our lab
 - The software has been discontinued and is no longer available for download
 - Does not support MAC OS

System Requirements

Required

2GHz CPU, 1GB RAM
50GB Hard Drive
Windows XP, Vista, Windows 7
Minimum 800x600 monitor

Recommended

3GHz CPU, 2GB RAM
50GB Hard Drive
Windows XP, Vista, Windows 7
Minimum 1280x1024 monitor

You can download a copy of NanoScope Analysis from
<ftp://Release:Z3Analysis91@sboftp.bruker-nano.com>. Once the software is installed, hit the F1 key or click on the blue button with the question mark on the toolbar for the online manual. A list of major new features is available under the "Help" menu – click on "What's New..."

No longer available for download




Propose solution

- The goal of this project is to develop an AFM image analysis program that can calculate the surface roughness average, Ra, using python
 - Ra is the arithmetic average of the absolute values of the profile heights over the evaluation length
 - Two modes: line Ra & area Ra
- AFM output data is a topography picture → extract the depth values from the picture
 - Calculate roughness average, Ra using eqn 1.

$$R_a = \frac{1}{L} \int_0^L |Z(x)| dx \text{ (eqn 1.)}$$


Where $Z(x)$ is the function that describes the surface profile analyzed in terms of height (Z) and position (x) of the sample over the evaluating length “ L ”






Overall workflow


User input





Read image
&
Convert to grayscale



Select region of interested



Output the Ra value and a height profile



Press the
“return” key on
keyboard



Method – Python modules

```
import cv2 # module import name for opencv
import numpy as np
from scipy import interpolate
import matplotlib.pyplot as plt
# generates an ipywidget that renders plots in a control
%matplotlib widget
from mpl_toolkits.mplot3d import Axes3D # 3D graphing
```

- cv2 → opencv
 - cv2.imread → read an image
- Numpy → mathematical operations on arrays
- Interpolate → linear interpolation
- Matplotlib → 2d plotting for heights profile
- Axes3D → 3d graphing for heights profile



Method – User inputs & image read

```
##### INPUT #####
imageName = 'pla3d0325.002_1'+'.png'


high = 1.5 # unit micron
low = 0 # unit micron

mode = 'line' #line or area

#####
# image reading
img = cv2.imread('/Users/huitingluo/CME502_Final_Project/AFM images/'+imageName)



# convert to Grayscale image
imgGray = cv2.cvtColor(img, cv2.COLOR_BGR2GRAY)

# ratio of scale bar and the range of grayscale values
ratio = (high-low)/(255-0)
```





Method – points selection



FAR
BEYOND

```
# Function: Pick Starting & ending points
def Pick_two_points(image, *args):
    global counter
    counter = 0
    # Create point matrix get coordinates of mouse click on image
    point_matrix = np.zeros((2,2),int)
    clone = image.copy()
    img = image.copy()
    def mousePoints(event,x,y,flags,params):
        global counter
        # Left button click
        if event == cv2.EVENT_LBUTTONDOWN: # indicates that the left mouse button is down
            if counter <2:
                point_matrix[counter] = x,y
                counter = counter + 1



        for x in range (0,2):
            #Draw circle on image with mouse click
            cv2.circle(img,(point_matrix[x][0],point_matrix[x][1]),8,(255,0,0),cv2.FILLED)
            cv2.putText(img,
                        'A' if x==0 else 'B',
                        (point_matrix[x][0]+10,point_matrix[x][1]-10),
                        cv2.FONT_HERSHEY_DUPLEX,
                        3.0,
                        (255,0,255),
                        3
                    )
        if counter == 2:
            # assign starting & ending points
            starting_x = point_matrix[0][0]
            starting_y = point_matrix[0][1]

            ending_x = point_matrix[1][0]
            ending_y = point_matrix[1][1]

            if 'line' in args:
                # Draw line between starting & ending points
                ImgWithMarker=cv2.line(img, (starting_x, starting_y), (ending_x, ending_y), (0, 255, 0), 3)
            if 'area' in args:
                # Draw rectangle for area of interest
                ImgWithMarker=cv2.rectangle(img, (starting_x, starting_y), (ending_x, ending_y), (0, 255, 0), 3)
```



Method – points selection



```
# keep the image window active |
while True:
    cv2.namedWindow('AFM', cv2.WINDOW_NORMAL)
    # Mouse click event on original image
    cv2.resizeWindow('AFM', image.shape[1], image.shape[0])
    cv2.setMouseCallback('AFM', mousePoints)
    cv2.imshow('AFM', img)
    key = cv2.waitKey(1) & 0xFF # maintain output window until user presses a key

    # if the 'r' key is pressed, <<<reset>>> the region
    if key == ord("r") and counter ==2:
        counter = 0
        point_matrix = np.zeros((2,2),int)
        img = clone.copy()
    # if the 'enter' key is pressed, break from the loop
    elif key == ord("\r"):
        break

    print('Done')
#    cv2.destroyAllWindows()
#    cv2.waitKey(1) ##### should destroy or not
#                                     #confirm selection
return point_matrix

#####
#####
```



Method – calculate line Ra

```
StartEndPoints = Pick_two_points(img, mode)

start_x = StartEndPoints[0][0]
start_y = StartEndPoints[0][1]
end_x = StartEndPoints[1][0]
end_y = StartEndPoints[1][1]

d_x = abs(end_x - start_x) # d_x : number of pixels at x axis
d_y = abs(end_y - start_y) # d_y : number of pixels at y axis

L = int(np.sqrt(d_x**2+d_y**2)+1) # L: Number of pixels in the line |

z_p = imgGray[min(start_y,end_y):max(start_y,end_y), min(start_x,end_x):max(start_x,end_x)]
z = high_ratio*(255-z_p)
# z = imgGray[min(start_y,end_y):max(start_y,end_y), min(start_x,end_x):max(start_x,end_x)] # start & end
#z = zg*100/255+50


if mode == 'line': # line interpolation

    f = interpolate.interp2d(np.arange(z.shape[1]), np.arange(z.shape[0]), z, kind='linear')
    #f = interpolate.interp2d(x, y, z, kind='linear')

    Xnew = np.linspace(0,z.shape[1],L)
    Ynew = np.linspace(0,z.shape[0],L)
    Znew = f(Xnew,Ynew)
    if (start_x-end_x)*(start_y-end_y)<0:
        Z_dial = np.fliplr(Znew).diagonal()
        if start_y>end_y:
            Z_dial = Z_dial[::-1]
    else:
        Z_dial = np.diagonal(Znew)
    if start_x > end_x and start_y>end_y:
        Z_dial = Z_dial[::-1]
    plt.plot(Z_dial)

Ra = np.sum(np.abs(Z_dial-np.mean(Z_dial)))/Z_dial.size
```

Linear interpolation






Method – calculate area Ra

```
elif mode == 'area':  
    x = np.arange(z.shape[1])  
    y = np.arange(z.shape[0])  
    X,Y = np.meshgrid(x,y)  
    fig = plt.figure()  
    ax = Axes3D(fig)  
    ax.plot_surface(X, Y, z, cmap=plt.get_cmap('twilight'))  
    ax.set_title("3d")  
    ax.set_xlabel("x")  
    ax.set_ylabel("y")  
    ax.set_zlabel("z")  
    #plt.savefig('3D.png') # save fig in the Main Folder  
    plt.show()  
  
Ra = np.sum(np.abs(z-np.mean(z)))/z.size  
  
print('Ra Roughness: ', Ra, 'um')
```




Result – line Ra



Done
Ra Roughness: 0.25251784593576704 um



Figure 1







Stony Brook University



Figure 1


Result – area Ra



Ra Roughness: 0.21828772297761592 um


Research: effect of coated substrate for cell growth – spincast coated PLA & 3d printed PLA

Spincast sample 1



$R_a = 0.0496 \text{ um}$


Spincast sample 2



$R_a = 0.0106 \text{ um}$


Flat

3d printed sample 1



$R_a = 0.2276 \text{ um}$

3d printed sample 2




$R_a = 0.4109 \text{ um}$


Rough

Research: effect of coated substrate for cell growth – spin cast coated PLA & 3d printed PLA



Rough



Flat



eGFP



Cells grow better in a flat surface

Conclusion

- The program of surface roughness analysis for AFM is developed
- Some advantage of the program
 - Beneficial for all group members in our lab
 - Both window and Mac users friendly
 - Free for everyone
 - Analyze the surface roughness of any region of interested
 - In nanoScope, image cropping is needed to get the area Ra of an interested area, and the cropping can only be square. Only horizontal line Ra is available

